



Product Datasheet

NINJ1 antibody (orb389333)


Description	NINJ1 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NINJ1. Please contact us for the exact immunogen sequence. The peptide is available as orb13907.
Target	NINJ1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	16 kDa
Uniprot ID	O70131 , P70617 , Q92982
NCBI	004139
Dilution Range	WB: 1:100-600, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P staining of rat kidney tissue using NINJ1 antibody.



Immunohistochemical staining of paraffin-embedded tissue using NINJ1 antibody.



WB analysis of mouse spinal cord (lane 1) using NINJ1 antibody.