

Product Datasheet

Ki67 antibody (orb389335)

Description

Ki67 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Ki67. Please contact us for the exact immunogen sequence. The peptide is available as orb13949.

Target

Ki67

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

358 kDa

Uniprot ID

[P46013](#), [E9PVX6](#), [Q5RJM0](#)

NCBI


[001139438](#), [001145966](#)

Entrez
[4288](#)
Dilution Range


IHC-P: 1:100-500, IF/ICC: 1:100-500

Expiration Date


12 months from date of receipt.



IHC-P image of guinea pig lung tissue us...



Immunohistochemical staining of paraffin...



IF analysis of guinea pig rectum tissue ...