



Product Datasheet

CD28 antibody (orb389337)


Description	CD28 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD28. Please contact us for the exact immunogen sequence. The peptide is available as orb13952.
Target	CD28
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	25 kDa
Uniprot ID	P31042, P10747
NCBI	001230006
Dilution Range	WB: 1:100-1000, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of guinea p...



IHC-P staining of rat thyroid tissue usi...