



Product Datasheet

MAP2K7 antibody (orb389339)


Description	MAP2K7 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MAP2K7. Please contact us for the exact immunogen sequence. The peptide is available as orb13962.
Target	MAP2K7
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	47 kDa
Uniprot ID	Q8CE90, Q4KSH7, O14733
NCBI	001284484
Dilution Range	WB: 1:100-800, IHC-P: 1:100-600
Expiration Date	12 months from date of receipt.



IHC-P staining of rat stomach tissue using...



WB analysis of rat heart (lane 1), mouse...



IHC-P staining of guinea pig colon tissue...