



Product Datasheet

TAL1 antibody (orb389344)


Description	TAL1 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TAL1. Please contact us for the exact immunogen sequence. The peptide is available as orb14045.
Target	TAL1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	34 kDa
Uniprot ID	P17542, D3Z921
NCBI	001274276
Dilution Range	WB: 1:100-800, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat lung...



Western blot analysis of rat liver (lane...)