



Product Datasheet

SPIN antibody (orb389348)


Description	SPIN antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SPIN. Please contact us for the exact immunogen sequence. The peptide is available as orb14185.
Target	SPIN
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	29 kDa
Uniprot ID	Q61142 , Q9Y657 , Q4V8J7
NCBI	006708
Dilution Range	WB: 1:100-800, IHC-P: 1:50-300
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of guinea p...



WB analysis of mouse heart (lane 1), mou...



Immunohistochemical staining of rat stom...