


Product Datasheet

COX IV antibody (Biotin) (orb389362)

Description	Goat polyclonal antibody against COX IV conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	COX4I1, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 17kDa band observed in Human Skeletal Muscle lysates (calculated MW of 19.6kDa according to NP_001852.1).
Protein Sequence	QGLASKWDYEKNE.
Clonality	Polyclonal
MW	19.6; 15.4
Entrez	1327
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...