


Product Datasheet

ITGA3 antibody (orb389369)

Description	Goat polyclonal antibody to ITGA3
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ITGA3 (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in Human and Mouse Heart lysates (calculated MW of 117kDa according to NP_002195.1).
Protein Sequence	KSQPSETERLTDD
Clonality	Polyclonal
MW	117; 119; 102
Entrez	3675
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human and Mouse...