


Product Datasheet

KRT19 antibody (Biotin) (orb389370)

Description	Goat polyclonal antibody against KRT19 conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	cytokeratin 19 (aa285-298), Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human Colorectal cancer lysates (calculated MW of 44.1kDa according to NP_002267.2).
Protein Sequence	HTEQLQMSRSEVTD.
Clonality	Polyclonal
MW	44.1
Entrez	3880
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Colorectal...