


Product Datasheet

PTN antibody (Biotin) (orb389377)

Description	Goat polyclonal antibody against PTN conjugated to Biotin.
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Pleiotrophin, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 18kDa band observed in Mouse fetal Brain lysates (calculated MW of 18.9kDa according to Mouse NP_032999.1).
Protein Sequence	KEGKKQEKMLD.
Clonality	Polyclonal
MW	18.9; 18.9; 18.9
Entrez	5764
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse fetal Brain...