


Product Datasheet

TPI1 antibody (Biotin) (orb389379)

Description	Goat polyclonal antibody against TPI1 conjugated to Biotin.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Triosephosphate isomerase, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 26kDa band observed in Human Liver lysates (calculated MW of 26.7kDa according to NP_000356.1).
Protein Sequence	LKPEFVDIINAKQ.
Clonality	Polyclonal
MW	26.7; 30.8; 18.0
Entrez	7167
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...