


Product Datasheet

12 Lipoxygenase/ALOX12 Antibody (orb389450)

Description	12 Lipoxygenase/ALOX12 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human 12 Lipoxygenase (186-231aa ALKRVTLLSSWNCLEDFDQIFWGQKSALAEKVRQCWQDDELFSYQ), different from the related mouse sequence by six amino acids, and from the related rat sequence by seven
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75694 MW
Uniprot ID	P18054
Dilution Range	Western blot, 0.1-0.5µg/ml, Mouse, Rat
Expiration Date	12 months from date of receipt.



WB analysis of
12
Lipoxygenase
using ant...