


Product Datasheet

GAPDH Antibody (orb389477)


Description

GAPDH Antibody


Species/Host	Rabbit
Reactivity	Gallus, Human, Monkey, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human GAPDH recombinant protein (Position: N136-E335). Human GAPDH shares 95% and 94.5% amino acid (aa) sequence identity with mouse and rat GAPDH, respectively.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36053 MW
Uniprot ID	P04406
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat
Expiration Date	12 months from date of receipt.



WB analysis of GAPDH using anti-GAPDH an...



IHC analysis of GAPDH using anti-GAPDH a...



IHC analysis of GAPDH using anti-GAPDH a...