

Product Datasheet

IFNGR1 Antibody (orb389482)

Description

IFNGR1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human IFNGR1 (443-484aa QELITVIKAPTSFGYDKPHVLVDLLVDDSGKESLIGYRPTED), different from the related mouse sequence by seventeen amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

54405 MW

Uniprot ID


P15260

Dilution Range


Immunohistochemistry (Frozen Section), 0.5-1µg/ml
Immunocytochemistry, 0.5-1µg/ml
Western blot, 0.1-0.5µg/ml
Flow Cytometry, 1-3µg/1x10⁶ cells

Expiration Date

12 months from date of receipt.



Flow
Cytometry
analysis of
A549 cells us...



WB analysis of
IFNGR1 using
anti-IFNGR1 ...