

Product Datasheet

pHNF4 antibody (FITC) (orb3895)

Description	pHNF4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HNF4 alpha around the phosphorylation site of Ser304
Target	phospho-HNF4 (Ser304)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41235
Dilution Range	Flow-Cyt=1ug/Test,IF=1:50-200
Expiration Date	12 months from date of receipt.