


Product Datasheet

Human Afamin ELISA Kit (orb389522)

Description	Human Afamin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L22-N599
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and urine.
MW	69069 MW
Uniprot ID	P43652
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100000, concentration in serum is around 25 ug/ml.
Expiration Date	Please enquire.



Human
Afamin
PicoKine
ELISA Kit
standard...