


Product Datasheet

Human Klotho beta/KLB ELISA Kit (orb389524)

Description	Human Klotho beta/KLB ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F53-L997
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	119808 MW
Uniprot ID	Q86Z14
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 100-213 ng/ml.
Expiration Date	Please enquire.



Human
Klotho beta
PicoKine
ELISA Kit
sta...