


Product Datasheet

Mouse DTK/TYRO3 ELISA Kit (orb389525)

Description	Mouse DTK/TYRO3 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A31-S418
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	96208 MW
Uniprot ID	P55144
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1 ng/ml.
Expiration Date	Please enquire.



Mouse
DTK/TYRO3
PicoKine
ELISA Kit
stand...