

Product Datasheet

Human PRDX2/Peroxiredoxin-2 ELISA Kit (orb389528)

Description Human PRDX2/Peroxiredoxin-2 ELISA Kit

Reactivity Human

Range 312 pg/ml - 20,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: A2-N198

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 21892 MW

Uniprot ID **P32119**

Sensitivity 10 pg/ml

Dilution Range Dilution ratio of 1:5, concentration in serum and plasma is around 13 ng/ml.

Expiration Date Please enquire.



Human
Peroxiredoxin
2 PicoKine
ELISA Kit...