


Product Datasheet

Human FOLR2/Beta Fr ELISA Kit (orb389529)

Description	Human FOLR2/Beta Fr ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T17-N230
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	29280 MW
Uniprot ID	P14207
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum is around 5.5 ng/ml.
Expiration Date	Please enquire.



Human
FOLR2
PicoKine
ELISA Kit
standard ...