


Product Datasheet

Human TNFRSF25/DR3 ELISA Kit (orb389530)

Description	Human TNFRSF25/DR3 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q25-F201
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	45385 MW
Uniprot ID	Q93038
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 312 pg/ml.
Expiration Date	Please enquire.



Human
TNFRSF25/DR3
PicoKine ELISA
Kit st...