


Product Datasheet

Mouse CD117/c-kit ELISA Kit (orb389531)

Description	Mouse CD117/c-kit ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S25-P527
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	109343 MW
Uniprot ID	P05532
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 3 ng/ml.
Expiration Date	Please enquire.



Mouse
CD117/c-kit
PicoKine
ELISA Kit
sta...