


Product Datasheet

Mouse KLK1/Kallikrein 1 ELISA Kit (orb389535)

Description	Mouse KLK1/Kallikrein 1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I25-D261
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	28775 MW
Uniprot ID	P15947
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 50-80 ng/ml.
Expiration Date	Please enquire.



Mouse KLK1
PicoKine
ELISA Kit
standard c...