


Product Datasheet

Mouse CA9/Carbonic anhydrase 9 ELISA Kit (orb389536)

Description	Mouse CA9/Carbonic anhydrase 9 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q32-D390
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	47265 MW
Uniprot ID	Q8VHB5
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Mouse CA9
PicoKine
ELISA Kit
standard
cu...