




Product Datasheet

PKA R1 alpha antibody (orb38954)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PRKAR1A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42982
Uniprot ID	P10644
NCBI	NP_997636.1 , NP_997637.1 , NP_002725.1 , NP_001263219.1 , NP_001263218.1 , NP_001265362.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of in CEM cell lin...



Western blot analysis of 293 cell lysate...