

Product Datasheet

Calbindin Antibody (orb389549)

Description	Calbindin Antibody
Species/Host	Mouse
Reactivity	Gallus, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Full-length recombinant human protein. HGNC name for this protein is CALB1.
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P05937
Entrez	793
Expiration Date	12 months from date of receipt.