



Product Datasheet

CLIC4 antibody (orb389649)


Description	CLIC4 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CLIC4. Please contact us for the exact immunogen sequence. The peptide is available as orb14326.
Target	CLIC4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	28 kDa
Uniprot ID	Q9Y696 , Q9Z0W7 , Q9QYB1
NCBI	039234
Dilution Range	WB: 1:200-1000, IHC-P: 1:100-800
Expiration Date	12 months from date of receipt.



IHC-P staining of guinea pig heart tissue...



WB analysis of rat heart (lane 1), rat s...



Immunohistochemical staining of rat kidn...