



Product Datasheet

CXCL5 antibody (orb389660)


Description	CXCL5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CXCL5. Please contact us for the exact immunogen sequence. The peptide is available as orb14363.
Target	CXCL5
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	12 kDa
Uniprot ID	P42830, P97885, P50228
NCBI	002985
Dilution Range	WB: 1:200-1500, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P staining of rat ovary tissue using...



WB analysis of rat stomach (lane 1), rat...



Immunohistochemical staining of paraffin...