

Product Datasheet

ompC antibody (orb389661)

Description	ompC antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from Escherichia coli IHE3034 OMPC (21-120/367 aa)
Target	OMPC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P06996
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.