

Product Datasheet

MPI-5a (orb389731)

Description	MPI_5a is a potent selective inhibitor of HDAC6, and poorly blocks other HDAC enzymes.
Note	For research use only
Formula	$C_{16}H_{17}N_3O_3$
Purity	> 98%,Standard References
MW	299.33
CAS Number	[1259296-46-2]
SMILES	<chem>CN1C=CC=C1C(=O)N2CCC3=C(C2)C=CC(=C3)C(=O)NO</chem>
Expiration Date	12 months from date of receipt.