

Product Datasheet

SQSTM1 antibody (orb390062)

Description Rabbit polyclonal antibody to SQSTM1.

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic peptide of human SQSTM1/p62(NP_003891.1).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 62kDa


Uniprot ID **Q13501**

Dilution Range WB: 1:500-2000


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat live...



Western blot analysis of extracts of 3T3...