



Product Datasheet

CEACAM5 antibody (orb390090)


Description	Rabbit polyclonal antibody to CEACAM5.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein of human CEACAM5
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	150kDa
Uniprot ID	P06731
Entrez	1048
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of mouse brain tissue.



Western blot analysis of extracts of mouse pancreas.