



Product Datasheet

HAPLN1 antibody (orb390093)


Description	Rabbit polyclonal antibody to HAPLN1.
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human HAPLN1(NP_001875.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	P10915
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue



Immunohistochemical staining of human esophageal tissue



Immunohistochemical staining of rat liver tissue