



Product Datasheet

HEXA antibody (orb390101)


Description	Rabbit polyclonal antibody to HEXA.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human HEXA(NP_000511.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19kDa/60kDa
Uniprot ID	P06865
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of mouse spleen tissue.



Immunohistochemical staining of mouse testis tissue.