



Product Datasheet

REG3A antibody (orb390117)


Description	Rabbit polyclonal antibody to REG3A.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human REG3A(NP_002571.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14-15kDa
Uniprot ID	Q06141
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of rat lung...



Immunohistochemical staining of mouse br...