



Product Datasheet

PI3 antibody (orb390118)


Description	Rabbit polyclonal antibody to PI3.
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human PI3(NP_002629.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12kDa
Uniprot ID	P19957
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of rat brain tissue.