



Product Datasheet

PRKAB1 antibody (orb390121)


Description	Rabbit polyclonal antibody to PRKAB1.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human AMPK β 1(NP_006244.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38kDa
Uniprot ID	Q9Y478
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunohistochemical staining of human es...



Western blot analysis of extracts from n...