



## Product Datasheet

### RPN1 antibody (orb390123)


|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RPN1.                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                      |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | IHC, WB                                                |
| <b>Immunogen</b>           | Recombinant fusion protein of human RPN1(NP_002941.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.     |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 68kDa                                                  |
| <b>Uniprot ID</b>          | <b>P04843</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human liver tissue.



Immunohistochemical staining of human tissue.