


Product Datasheet

TCEB1 antibody (orb390131)

Description	Rabbit polyclonal antibody to TCEB1.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TCEB1(NP_001191791.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12kDa
Uniprot ID	Q15369
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...