


Product Datasheet

SUMO1 antibody (orb390136)

Description	Rabbit polyclonal antibody to SUMO1.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide of human SUMO1(NP_003343.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8kDa/11kDa
Uniprot ID	P63165
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...