



Product Datasheet

WISP2 antibody (orb390141)


Description	Rabbit polyclonal antibody to WISP2.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human WISP2(NP_003872.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa/26kDa
Uniprot ID	076076
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of human tissue.



Immunohistochemical staining of human tissue.