



Product Datasheet

KLK11 antibody (orb390148)


Description	Rabbit polyclonal antibody to KLK11.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human KLK11(NP_001129504.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	Q9UBX7
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of rat kidney tissue.