



Product Datasheet

CBX5 antibody (orb390159)


Description	Rabbit polyclonal antibody to CBX5.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human HP1 alpha/HP1 alpha/CBX5(NP_001120794.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa
Uniprot ID	P45973
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of human breast tissue.



Immunohistochemical staining of human breast tissue.