



Product Datasheet

EED antibody (orb390201)


Description	Rabbit polyclonal antibody to EED.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human EED(NP_003788.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50-75kDa
Uniprot ID	075530
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunofluorescence analysis of L929 cells...



Western blot analysis of extracts of various cell lines...