


Product Datasheet

NOS1 antibody (orb390203)

Description	Rabbit polyclonal antibody to NOS1.
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human NOS1(NP_000611.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150kDa
Uniprot ID	P29475
Dilution Range	WB: 1:500-2000, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Rat...