


Product Datasheet

RCOR1 antibody (orb390221)

Description	Rabbit polyclonal antibody to RCOR1.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human CoREST/CoREST/RCOR1(NP_055971.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa
Uniprot ID	Q9UKL0
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...