



Product Datasheet

FN1 antibody (orb390247)


Description	Rabbit polyclonal antibody to FN1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human Fibronectin (NP_002017.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	272KDa
Uniprot ID	P02751
Entrez	2335
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse le...



Immunohistochemical staining of human li...



Immunohistochemical staining of human lu...