

Product Datasheet

S1PR1 antibody (orb390249)

Description	Rabbit polyclonal antibody to S1PR1.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human EDG1/HEXIM1(NP_001391.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	P21453
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.

