



Product Datasheet

ECH1 antibody (orb390255)


Description	Rabbit polyclonal antibody to ECH1.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ECH1(NP_001389.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa
Uniprot ID	Q13011
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human liver tissue.



Immunohistochemical staining of human colon tissue.