


## Product Datasheet

### NR1I2 antibody (orb390271)


|                            |                                                            |
|----------------------------|------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NR1I2.                       |
| <b>Species/Host</b>        | Rabbit                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                          |
| <b>Conjugation</b>         | Unconjugated                                               |
| <b>Tested Applications</b> | WB                                                         |
| <b>Immunogen</b>           | A synthetic peptide of human PXR (NP_071285.1).            |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                |
| <b>Note</b>                | For research use only                                      |
| <b>Isotype</b>             | IgG                                                        |
| <b>Clonality</b>           | Polyclonal                                                 |
| <b>Purity</b>              | Affinity purification                                      |
| <b>MW</b>                  | 50KDa                                                      |
| <b>Uniprot ID</b>          | <b>075469</b>                                              |
| <b>Entrez</b>              | <b>8856</b>                                                |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                            |



Immunohistochemical staining of rat panc...



Immunohistochemical staining of rat test...



Western blot analysis of extracts of rat...